

I-Bridge®							
Implant system	Platform	I-Bridge®	I-Bridge® Zirconia	I-Bridge® 2	I-Bridge® Evolution	I-Pack® available	I-Flex™ available
BRÄNEMARK SYSTEM® Implant Level	NP	●		●	●	914510	914515
	RP	●		●	●	914110	914115
	WP	●	●			914410	
BRÄNEMARK SYSTEM® Multi Unit Abutment Level	NP	●	●	●	●	914100	*
	RP	●	●	●	●	914100	*
	WP	●	●			914400	
Straumann® Tissue Level	NN						
	RN	●		●	●	914650	914655
	WN	●		●		914660	*
Straumann® SynOcta® Abutment Level	NN						
	RN	●				914620	
	WN						
Straumann® Bone Level	NC 3.3	●		●		914670	911015
	RC 4.1	●		●		914680	911025
	RC 4.8	●		●		914690	911035
NobelReplace® Implant Level	NP	●		●		914520	914525
	RP	●		●		914120	914125
	WP	●		●		914420	*
NobelReplace® Multi Unit Abutment Level	NP	●	●	●	●	914100	*
	RP	●	●	●	●	914100	*
	WP	●	●	●	●	914100	*
Astra Tech Implant System® Implant Level	3.5/4.0	●		●		914850	914815
	4.5/5.0	●		●		914870	914825
Astra Tech Implant System® Abutment Level	20°	●		●		914810	*
	45°	●		●		914830	*
Neoss™	3.5	●		●		916120	914015
	4.0	●		●		916120	914015
	4.5	●		●		916120	914015
3i Certain®	3.4	●		●		916320	913015
	4.1	●		●		916330	913025
	5.0	●		●		916340	913035
BioHorizons Internal	3.5	●		●		916420	911115
	4.5	●		●		916430	911125
	5.7	●		●		916440	911135

I-Bridge®
Implant system

	Platform	I-Bridge®	I-Bridge® Zirconia	I-Bridge® 2	I-Bridge® Evolution	I-Pack® available	I-Flex™ available
BioHorizons External	3.5	●		●		916450	911145
	4.0	●		●		916460	911155
	5.0	●		●		916470	911165
	6.0	●		●		916480	911175
Zimmer Dental Tapered Screw-Vent®	3.5	●		●		916420	911115
	4.5	●		●		916430	911125
	5.7	●		●		916440	911135
Camlog® Implant System	3.3	●		●		914940	*
	3.8	●		●		914941	*
	4.3	●		●		914942	*
	5.0	●		●		914943	*
NobelActive™	3.5	●		●			*
	4.3	●		●			*
Ankylos® Balance C Abutment Level	4.2	●		●			*
	5.5	●		●			*

* Mark screw hole manually by using an Angulation Guide.

Version: 2011-04-18

BRÄNEMARK SYSTEM® and NobelReplace® are registered trademarks of Nobel Biocare AB, Sweden. Straumann® is a registered trademark of Straumann Holding AG, Schweiz.

Astra Tech Implant System® is a registered trademark of Astra Tech AB, Sweden. Neoss™ is a trademark application of Neoss Ltd, United Kingdom. Certain® is a registered trademark of Biomet 3i Inc, USA.

Screw-Vent® is a registered brand of Zimmer Dental, USA. Camlog® is a registered trademark of CAMLOG Biotechnologies AG, Germany. NobelActive™ is a trademark application of Nobel Biocare AB, Sweden.

Ankylos® is a registered trademark of DeguDent GmbH, Germany. BioHorizons® is a registered trademark of BioHorizons Implant Systems Inc, USA.