

Heraeus Kulzer

Mitsui Chemicals Group

cara 

cara I-Bridge[®]

Overview of platforms and systems



cara makes life so easy



cara I-Bridge®

Implant System	Platform	cara I-Bridge regular			cara I-Bridge angled				
		Digital Design	CoCr milled	CoCr SLM	Ti	CoCr milled	CoCr SLM	Ti	cara I-Pack
Ankylos® Balance Base C Abutment system	4.2	—	●	●	●	●	●	●	*
	5.5	—	●	●	●	●	●	●	*
Astra Tech OsseoSpeed® Abutment Level	20°	—	●	●	●	●	●	●	*
	45°	—	●	●	●	●	●	●	*
Astra Tech OsseoSpeed® Switch Level	3.5/4.0	—	●	●	●	●	●	●	●
	4.5/5.0	—	●	●	●	●	●	●	●
BioHorizons® External Implant Level	3.5	—	●	●	●	●	●	●	●
	4.0	—	●	●	●	●	●	●	●
	5.0	—	●	●	●	●	●	●	●
	6.0	—	●	●	●	●	●	●	●
BioHorizons® Internal Implant Level	3.0	—	●	●	●	●	●	●	*
	3.5	✓	●	●	●	●	●	●	●
	4.5	✓	●	●	●	●	●	●	●
	5.7	✓	●	●	●	●	●	●	●
Biomet 3i™ Certain® Implant Level	3.4	✓	●	●	●	●	●	●	●
	4.1	✓	●	●	●	●	●	●	●
	5.0	✓	●	●	●	●	●	●	●
Biomet 3i™ External Implant Level	3.4	✓	●	●	●	●	●	●	●
	4.1	✓	●	●	●	●	●	●	●
	5.0	✓	●	●	●	●	●	●	●
	6.0	✓	●	●	●	●	●	—	*
Biomet 3i™ Abutment level	Low Profile Abutment	—	—	—	—	●	●	●	—
Brånemark System® Implant Level	NP	✓	●	●	●	●	●	●	●
	RP	✓	●	●	●	●	●	●	●
	WP	✓	●	●	●	—	—	—	●
Brånemark System® Multi Unit Abutment Level	NP	✓**	●	●	●	●	●	●	*
	RP	✓**	●	●	●	●	●	●	*
	WP	✓**	●	●	●	—	—	—	●



Implant System	Platform	cara I-Bridge regular			cara I-Bridge angled				
		Digital Design	CoCr milled	CoCr SLM	Ti	CoCr milled	CoCr SLM	Ti	cara I-Pack
Camlog® Implant System Implant Level	3.3	✓**	●	●	●	●	●	●	●
	3.8	✓**	●	●	●	●	●	●	●
	4.3	✓**	●	●	●	●	●	●	●
	5.0	✓**	●	●	●	●	●	●	●
Dentaurum tiologic® Implant Level	S	✓**	●	●	●	●	●	—	*
	M	✓**	●	●	●	●	●	—	*
	L	✓**	●	●	●	●	●	—	*
Dentaurum tiologic® Abutment Level	S	✓**	●	●	●	●	●	—	*
	M	✓**	●	●	●	●	●	—	*
	L	✓**	●	●	●	●	●	—	*
Neoss™ Implant Level	3.5	✓**	●	●	●	●	●	●	●
	4.0	✓**	●	●	●	●	●	●	●
	4.5	✓**	●	●	●	●	●	●	●
NobelActive™ Implant Level	3.5	✓	●	●	●	●	●	●	●
	4.3	✓	●	●	●	●	●	●	●
NobelReplace® Implant Level	NP	✓	●	●	●	●	●	●	●
	RP	✓	●	●	●	●	●	●	●
	WP	✓	●	●	●	●	●	●	●
NobelReplace® Multi Unit Abutment Level	NP	✓**	●	●	●	●	●	●	*
	RP	✓**	●	●	●	●	●	●	*
	WP	✓**	●	●	●	●	●	●	*
Osstem® GS/TS Abutment Level	Mini	—	●	●	●	●	●	●	*
	Standard	—	●	●	●	●	●	●	*
Osstem® SS Implant Level	Mini	—	●	●	●	●	●	●	*
	Standard	—	●	●	●	●	●	●	*
	Wide	—	●	●	●	●	●	—	*
Osstem® US Implant Level	Mini	—	●	●	●	●	●	●	*
	Standard	—	●	●	●	●	●	●	*
	Wide	—	●	●	●	●	●	●	*



cara I-Bridge®

Implant System	Platform	cara I-Bridge regular			cara I-Bridge angled				
		Digital Design	CoCr milled	CoCr SLM	Ti	CoCr milled	CoCr SLM	Ti	cara I-Pack
Straumann® Bone Level Implant Level	NC 3.3	✓	●	●	●	●	●	●	●
	RC 4.1	✓	●	●	●	●	●	●	●
	RC 4.8	✓	●	●	●	●	●	●	●
Straumann® Bone Level Multi-Base Abutment Level	3.5	—	●	●	●	●	●	—	*
	4.5	—	●	●	●	●	●	—	*
Straumann® Implant Level	NN	—	—	—	—	●	●	●	—
	RN	✓	●	●	●	●	●	●	●
	WN	✓	●	●	●	●	●	●	●
Straumann® SynOcta Abutment Level	NN	—	—	—	—	—	—	—	—
	RN	—	●	●	●	●	●	●	—
	WN	—	—	—	—	—	—	—	—
Thommen SPI® Implant Level	3.5	✓**	●	●	●	●	●	●	●
	4.0	✓**	●	●	●	●	●	●	●
	4.5	✓**	●	●	●	●	●	●	●
	5.0	✓**	●	●	●	●	●	●	●
	6.0	✓**	●	●	●	●	●	●	●
Thommen® Abutment level	VARIOMulti	—	●	●	●	●	●	●	—
TRI®-Vent Implant Level	3.5	—	●	●	●	●	●	—	*
Xive®/Frialit® Implant Level	3.4	✓	●	●	●	●	●	●	*
	3.8	✓	●	●	●	●	●	●	*
	4.5	✓	●	●	●	●	●	●	*
	5.5	✓	●	●	●	●	●	●	*
Zimmer Dental Tapered Screw-Vent® comp.	3.5	✓	●	●	●	●	●	—	*
	4.5	✓	●	●	●	●	●	—	*
	5.7	✓	●	●	●	●	●	—	*

cara Information

* The screw channel's exit hole is marked with a raised indication, e.g. using wax. cara I-Flex with an Access guide or Angulation Guide shall be used for all analogues so as not to exceed the maximum angulation of the screw holes.

** cara I-Bridge can be designed virtually with other scanbodies than NT Trading or Medentica. Please contact Technical Hotline for further information.

cara I-Bridge new connection

We can now manufacture cara I-Bridge to any implant system. Please contact Customer Service for more information about how to order a cara I-Bridge new connection. A great opportunity that help customer to get a great prosthetic solution to all patient situations.

cara I-Bridge Zr

cara I-Bridge is available in Zr to the following platforms:

Regular/Angled:

- Brånemark system Implant level NP/RP
- Brånemark system Multi Unit Abutment level NP/RP
- Nobel Replace Multi Unit Abutment level NP/RP/WP
- Straumann system Implant level RN/WN

Regular:

- Brånemark system Implant level WP
- Brånemark system Multi Unit Abutment level WP
- Straumann system Implant level RN/WN

All products are not available in all markets.

Brånemark System®, NobelReplace® and NobelActive® are trademarks of the Nobel Biocare group. Straumann® are trademarks of Straumann Group. Astra Tech OsseoSpeed®, Xive®, Frialit® are trademarks of Dentsply Friadent Group. 3i®, Biomet®, Biomet 3i™, Certain® are trademarks of Biomet 3i, Inc. Screw-Vent® is a registered brand of Zimmer Dental, USA. Thommen® Implant System are trademarks of Thommen Medical Group.

All trademarks and company names are the property of their respective owners.

Digital Design

The customer designs cara I-Bridge digitally, using their own experience and adapting to the patient's specific requirements virtually. The file is then uploaded and sent for production. This flexibility reduces time and cost for both laboratories and clinicians.

The customer can use either 3Shape, sending the file via Heraeus Kulzers' inbox, or via open STL-files uploaded to Meeting Point.

The leading scanners/software in the market; 3Shape, Exocad and Dental Wings can be used for the design.

For more information, please visit www.cara-kulzer.com

Design recommendation for 3Shape

The customer will need the following to be able to start

- 3Shape DentalManager 2013, Ver. 2.8.8.8.
- Updated library for the abutments in 3Shape from Heraeus Kulzer.
- Scan body; We support only scanbodies from Medentica and NT Trading in our material file. If other scanbodies are used please contact Technical Hotline for information.



Our specialists in the field will be pleased to provide you with information and service.

ORDER & CUSTOMER SERVICE

☎ +46 42 453 07 00

TECHNICAL HOTLINE

☎ +46 42 453 07 00

GIVE US A CALL!

www.cara-kulzer.com



List of our distributors upon request



Contact in Nordic & Baltic States

Heraeus Kulzer Nordic AB
Box 22237
SE-250 24 Helsingborg
Tel: +46 42 453 07 00
Fax: +46 42 453 07 01
cara-service@kulzer-dental.com
www.cara-kulzer.se



Contact in Germany

Heraeus Kulzer GmbH
Grüner Weg 11
63450 Hanau
Tel: 0800 437 25 22
Fax: 0800 437 23 29
cadcam@kulzer-dental.com
www.cara-kulzer.de

Contact in UK

Heraeus Kulzer Limited
Heraeus House, Albert Road
Newbury, Berkshire
RG14 1DL
www.cara-kulzer.com



Intertek